



# Biogel<sup>®</sup> M


*Formulated for improved grip*



- Formulation identical to Biogel surgical glove, but with micro-textured surface for enhanced grip
- Unique shape for hard-to-fit hands
- Powder-free to eliminate starch powder-related complications
- Biogel coating for ease of donning with damp or dry hands
- 1.0 AQL for holes<sup>1</sup>
- 100% air inflation and visually inspected<sup>2</sup>
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging

## Product Properties

- Material: Natural Rubber Latex with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Nominal thickness: fingertip 0.26mm (10mils)
- Straight finger design
- Non-pyrogenic
- Low dermatitis potential<sup>3</sup>
- Contains 50 µg or less of water extractable protein per gram<sup>4</sup>
- Color: Straw



Order #	Size	Pairs/Box
30555	5½	50
30560	6	50
30565	6½	50
30570	7	50
30575	7½	50
30580	8	50
30585	8½	50
30590	9	40

4 Boxes per Case