



# Biogel® Skinsense®


*A non-latex surgical glove made of chloroprene*



- Chloroprene composition eliminates possibility of glove-related latex protein sensitization
- Powder-free to eliminate starch powder-related complications
- Biogel coating for ease of donning with damp or dry hands
- 1.0 AQL for holes<sup>1</sup>
- 100% air inflation and visually inspected<sup>2</sup>
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging

## Product Properties

- Material: Chloroprene with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Non-pyrogenic
- Low dermatitis potential<sup>3</sup>
- Nominal thickness: fingertip 0.20mm (8mils)
- Curved finger design
- Color: Neutral tone



Order #	Size	Pairs/Box
31455	5½	50
31460	6	50
31465	6½	50
31470	7	50
31475	7½	50
31480	8	50
31485	8½	50
31490	9	40

4 Boxes per Case

<sup>1</sup>1.0 AQL for freedom from holes ASTM 3577, Test Method D5151.

<sup>2</sup>Each Biogel glove undergoes multiple checks for quality and user confidence.

<sup>3</sup>Low dermatitis potential – gloves demonstrated reduced potential for sensitizing user to chemical additives.