



Biogel®


The World's Finest Surgical Gloves.™



- Acknowledged by many as the finest general-use surgical glove
- Micro-roughened surface
- Outstanding fit, feel and comfort
- Powder-free to eliminate starch powder-related complications
- Biogel coating for ease of donning with damp or dry hands
- 1.0 AQL for holes¹
- 100% air inflation and visually inspected²
- Meet or exceed ASTM surgical glove standards
- Beaded cuff
- Sure-Flap packaging

Product Properties

- Material: Natural Rubber Latex with Biogel polymer inner coating
- Biogel coating soothes the skin to help prevent moisture loss
- Nominal thickness: fingertip 0.26mm (10mils)
- Curved finger design
- Non-pyrogenic
- Low dermatitis potential³
- Contains 50 µg or less of water extractable protein per gram⁴
- Color: Straw



Order #	Size	Pairs/Box
30455	5½	50
30460	6	50
30465	6½	50
30470	7	50
30475	7½	50
30480	8	50
30485	8½	50
30490	9	40

4 Boxes per Case

¹ 1.0 AQL for freedom from holes ASTM 3577, Test Method D5151.

² Each Biogel glove undergoes multiple checks for quality and user confidence.

³ Low dermatitis potential – gloves demonstrated reduced potential for sensitizing user to chemical additives. **WARNING: Do not use this product if you have a known allergy to natural rubber protein or chemical additives.**

⁴ This product contains 50 µg or less of water extractable protein per gram. **CAUTION: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.**

