

Monopolar L-hook electrode

Truly customised laparoscopic solutions



The 5mm Monopolar L-hook electrode enables the precise, controlled dissection and coagulation of vessels and tissues. Available both as sterile single-packed product and for inclusion in Mölnlycke® procedure trays.

- Convenient single-use instrument guarantees sterility and effective monopolar current containment, preventing adjacent site tissue burn
- Lightweight, ergonomic handle design allows precise controlled tissue dissection
- Fully-insulated 33cm shaft provides safe current containment
- Available with single-patient-use cable, for optimal flexibility and convenience



Monopolar L-hook electrode without cable



Hook



Monopolar L-hook electrode with cable

Ref. No.	Description	Size	Pouch size	Pcs/disp.
899100	Monopolar L-hook electrode	5mm	33cm	10
899101	Monopolar L-hook electrode with cable	5mm	33cm	10

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00
 The Mölnlycke trademark, name and logo is registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 © 2017 Mölnlycke Health Care AB. All rights reserved. HQSU0324